



Products:

Electronic Lab Notebook (ELN) - Designed to replace paper notebooks for Research, Analytical and Production laboratories.

Laboratory Information Management (LIMS) – and easily implemented and cost effective solution to sample & results tracking.

Project Management Dashboard (PMD) - a collaborative Project Management tracking system that your whole organization can use.

Editions:

Team Editions are designed and priced attractively for smaller groups numbering from 5 to 200 users.

Enterprise Edition are designed for larger scale implementations numbering from 50 to 10,000 users. The Enterprise Version is developed to meet typical WAN bandwidth and global requirements.

Personal Editions install on a Windows PC allowing users to edit ELN notebooks locally. When combined with a Team or Enterprise version, users may check-in and check-out notebooks from the master database.

*Minimum System Requirements **

LABTrack Personal and Personal Lite Editions

Windows 8 Professional, Business or Ultimate or newer
(Home versions of Windows are not supported)
DuoCore Processor, 2.0 GHZ CPU or better
2 GB RAM, 20 GB free disk space
.NET Framework 3.5
IIS 7 and newer
Internet Explorer 9 and newer
SQL Express 2012 and newer

LABTrack Team Edition

Windows 2008 Server OS or newer
DuoCore Processor, 2.0 GHZ CPU or better
8 GB RAM or more, 20 GB free disk space
Microsoft SQL Server or SQL Express 2012 or newer
IIS 7 or newer
.NET Framework 3.5
AJAX Extensions (LIMS)
Internet Explorer 9 and newer
Microsoft Reporting Services required for custom reports.
Application and Database Server can be hosted on the same server.
Additional disk space as required by database storage.

LABTrack Enterprise Edition

Application Server -

Windows 2008 Server OS or newer
DuoCore Processor, 2.0 GHZ CPU or better
8 GB RAM or more, 20 GB free disk space
Microsoft SQL Server or SQL Express 2012 or newer
IIS 7 or newer
.NET Framework 3.5
Internet Explorer 9 and newer
Microsoft Reporting Services required for custom reports.

Database Server – Windows 2008 Server OS or newer
and either SQL Server or SQL Express 2012 or newer

Additional disk space as required by database storage.

Team & Enterprise Client PC Requirements

Windows 8 or newer or Mac X OS or newer
2 GB RAM, 50 MB free disk space
Internet Explorer, Firefox, Chrome or Safari Browser

* Requirements are subject to change without notice.

LABTrack, LLC

22365 El Toro Rd #185, Lake Forest, CA 92630
Tele: (US) +1 (949) 581-9004, www.labtrack.com